



**AUT-701** 




AUT-701
Unique & Powerful
Titration System

Powerful functions can strongly support QA and R&D





Adding a second burette and other items to the standard configuration (single titration system) provides capability to perform two different titration methods (such as neutralization and redox) at the same time. Simultaneous measurement of pH and neutralization titration is also possible. By connecting two Model TTT-510 turntables, simultaneous measurement of multiple samples is also possible.




## Example of dual titration system (pH & ORP)

Standard (single) titration system with additional equipment added (ABT-7 Auto Burette, ST-7 Magnetic Stirrer, ORP Electrode, Electrode Stand & Electrode Holder).



## Back lit Color LCD provides Excellent Ledgibility



#### Complete Range of Titration Methods Available (up to 4 electrodes/ 10 burettes)

Standard configuration includes 2pH/mV inputs. Optional titration units (photometric, polarization, conductivity, potentiometric) can be connected to the third and forth units. In addition, up to 10 sets of burette can be connected. Multi-titration and advanced sequential mode titrations are also available.

#### Advanced validation support function (display and logging of ambient temperature & RH)

pH calibration history (max 20 calibrations) function, burette volume checking history function and logging of ambient temperature and relative humidity using optional temp/RH sensor can offer advanced quality control.





Temp/RH sensor

# USB & LAN Interfaces Available as standard

Data acquisition software, 701, can collect and store data in CSV format for easy data management and transfer to other software applications.

## **Enriched Data Memory Function**

100 data sets can be re-analyzed for single titration system.

(200 data sets for dual system)

#### **Fast Titration Mode**

Titration time can be shortened. (Not applicable for some applications.)

# Electrode & Titration Kit can be Applied to Various Titrations

#### Micro-titration kit (option)

Micro-titration kit is ready with combined pH electrode for micro solution and silver combined electrode.


Waste water can be reduced.

(Not applicable for some applications.)

Minimum volume: approx. 10 mL Maximum volume: approx. 20 mL

#### ●Content

Adaptor, Titration nozzle,
Titration cell (50 pcs.),
O-ring P-5 (2 pcs.), O-ring
P12 (2 pcs.) ,Stirring bar (2 pcs.)



### Worry-free operation with durable glass electrode


Strong & Float-pH electrodes

Strong: Tip part is reinforced.
(more than ten times tougher than our


Float : Built-in float enables recognition of inner solution concentration

conventional type)


at a glance.



## Print sample



# Modular Concept Allows Expansion for Your Future Needs




# System Examples

Automatic Multiple Sample Titration System

#### AUT-701 + TTT-510


Multiple sample titration can be done automatically by setting pretreated samples on the Turntable. Flex function allows operator to change titration condition of each sample. (For example, Table #1 to #5: pH measurement, over #6: neutralization titration)



# Dual Automatic Multiple Sample Titration System

#### AUT-701(Dual System) + TTT-510 × 2


Different two kinds of multiple titration system can be configured.



# Automatic Multiple Sample Titration System with Automatic Sampler

#### AUT-701 + TTT-510 + AST-3210

Series of multiple sample operations (measurement - dilution - reagent dispense- titration - cleaning) can be performed automatically after setting samples on the Turntable.



# Peripheral Equipment



## **Automated Continuous Multiple Sample Measurement Capability**

# Turntable TTT-510

- Chemical resistant cleaning reservoir
- ●Various electrode cleaning methods
- \*\*Optional air pump, reagent pump and others are required for above cleaning methods except for pure water shower.

#### Electrode holding tank available for overnight operation

After the final measurement, the electrode is stored in a storage tank filled with pure water, so there are no concerns that the electrode will dry out.

#### **Specifications**

Display unit	LCD display		
Operating keys/switch	Flat keyboard		
Number of samples /beakers used	12 samples 200 mL, 300mL tall beaker 18 samples 100 mL tall beaker (Manufacturer's designation) 36 samples 30 mL, 50 mL tall beaker (Manufacture's designation) 60 samples/(100 samples) 20 mL cell (Manufacturer's designation)		
Cleaning mode	① Pure water shower ② Pure water bubbling ③ Liquid chemical shower ④ Pure water shower ⑤ Liquid chemical bubbling → Pure water shower ⑤ Liquid chemical bubbling → Pure water shower → Air blow ⑥ Liquid chemical bubbling → Liquid chemical bubbling → Liquid chemical shower → Air blow ⑥ Liquid chemical shower → Air blow		
Cleaning time	1 – 9,999 seconds (Setting possible)		
Air blow frequency	1 – 9 times (Setting possible)		
End detection	Using endpin or setting final sample number by key operation		
Cleaning tank material	Teflon, Daiflon		
Electrode storage tank	Standard		
Stirring time before measurement	0 - 9,999 seconds		
Waiting time before measurement	0 - 9,999 seconds		
Alarm	Cleaning water 1 (Pure water) Empty, Cleaning water 2 (Chemical) Empty, Waste water Full		
Ambient temperature	5 − 40℃		
Power source	AC100V - 240V 50/60Hz		
Power consumption	Max. approx. 55VA (without optional pump & valve) Max. approx. 175VA (with optional pump & valve)		
External dimensions	Main unit : Approx. W380 mm × H392 mm × D622 mm  Maximum size when moving : Approx. W 415mm (with table plate) × H518 mm × D622 mm		

#### Standard Accessories

Table (designated for 12, 18, 36 and 60 samples)		
Electrode holder (designated samples / use)		
Pure water tank (10 L, with float switch)		
Waste water tank (10 L, with float switch)		
Pure water tube		
Waste water tube		
End pin		
Polyethylene Beaker (200 mL)		
Instruction Manual		

 $<sup>\*</sup>$  The attachment cable for the titrator is sold separately based on the model designation.

#### **Options**

Description	Code
Connection cable (2m)	7075670K
Air pump	P000023
Liquid chemical pump	P000022
Bubbling pump	OFA00001
Waste water valve (pinch valve)	OQA00001
Waste water valve (solenoid valve)	0QA00002

#### Full Automatic Operation ~ Measurement to Cleaning

# Automatic Sampler AST-3210



Sample measuring principle	Three-way stopvalve switching measuring tube
Sample volume	About 0.2 – 20 mL (fixed loop volume [measuring tube])
Measurement repeatability	0.5% or less for 10 mL sample as CV value
Measurement vessel	Material : hard glass, Number of electrode : Max. 4 Max. volume : approx.150 mL, Number of titration nozzle : Max. 4 Min. required volume : about 100 mL
Alarm indication	Pure water empty, Waste water full, System abnormal
Tanks for waste water & pure water	20 L plastic tank
Power source	AC line 50/60Hz
Power consumption	Approx. 100VA
External dimensions & Weight	Approx. W565 mm $\times$ H635 mm $\times$ D480 mm, approx. 45 kg

<sup>\*</sup>The attachment cable for the titrator, trunk cable for electrode and others are sold separately based on the model designation

#### For System Expansion (For Example Dual Titration)



Magnetic Stirrer ST-7



#### **Data Acquisition Software for PC**

# Data Acquisition Software 701 7075650K

Capable of importing data to PC in the CSV format using USB memory, LAN or RS-232.

#### For Plain Paper Printing

External printer EPS-R



# **Optional Units**

# Photometric Titration Unit Two wave length auto switching available

FUT-7010 (Immersion electrode OPE-21A attached)

Used to detect the color variation of indicator and used for quantitative analysis of metallization like plating solution.

\*The adaptor is built into the titrator and the amp. unit is standalone.



Amp. unit dimension : Approx. W102 X H153 X D119 mm Standard interference filter : 530 nm, 630 nm

#### **Polarization Titration Unit**

FUT-7030 (Twin platinum electrode HDP-303 attached)

Used for iodide value or diazo titration.

\*This unit is built into the titrator.

Measuring method : constant-current voltage method, constant-voltage current method

Detecting electrode: twin platinum electrode (HDP-303) Applied voltage: 0 +/- 1,000 mV (settable at optimum voltage) Applied current :  $0 + /- 1,000 \mu A$  (settable at optimum voltage)



#### **Potentiometric Titration Unit**

FUT-7050 (Electrode is option)

\*This unit is built into the titrator.

mV : 0 +/- 2,000.0 mV Temperature : 0-100.0℃



#### **Electric Conductivity Titration Unit**

FUT-7040 (Conductivity cell CT-57101B attached)

Can be used for electric conductivity measurement.

\*This unit is built into the titrator.

Measuring method: AC two electrode method

Temperature

: 0-200.0S/m Measuring range (depends on cell)

Manual ranging when titration set by titrator.

Auto ranging when electric conductivity measurement

: 0-100.0℃

#### Electrodes

Application	Description	Code	Remarks
	pH combined electrode (general purpose)	GST-5741C	One of standard accessories of AUT-701
Neutralization titration	pH combined electrode (precision measurement)	ELP-035	
	pH combined electrode (double junction)	ELP-062	Junction is replaceable. Effective against KCI inner solution outflow problem
Redox titration	ORP combined electrode (general purpose)	PS-5111C	
Redox (Itration	ORP combined electrode (double junction)	ELM-016	Junction is replaceable. Effective against KCI inner solution outflow problem
Salinity titration	Silver combined electrode (general purpose)	ELX-006	Junction is replaceable
Non-aqueous titration	Glass electrode (general purpose)	HGS-2005	Optional electrode adapter (OJD00001) and electrode holder
Non-aqueous titration	Reference electrode (double junction)	HS-305DS	(S-HLD-S) will be required.
Photometric titration	Immersion electrode	OPE-21A	Attached with photometric titration unit
Polarization titration	Twin platinum electrode	HPD-303	Attached with polarization titration unit
Electric conductivity titration	Electric conductivity cell	CT-57101B	Attached with electric conductivity titration unit

#### Parts & Solutions

#### Parts for titrator

Description	Code	Remarks
Magnetic stirrer	ST-7	
Electrode stand with stopper	6948810K	
Electrode holder	OIB00001	Standard accessories of AUT-701
Electrode attachment (G)	0IB00004	For dual system separate order would be required
Electrode attachment (J)	0IB00005	
Electrode attachment (N)	0IB00008	
Electrode adapter	OJD00001	Required for using glass electrodes and reference electrodes
Electrode holder	S-HLD-S	hequired for using glass electrodes and reference electrodes
Electrode relaying lead wire (1m)	OGB00001	
Electrode relaying lead wire (3m)	0GB00002	
Micro-titration kit	7075600K	Content: Adaptor, Titration nozzle, Titration cell (50 pcs.), O-ring P-5 (2 pcs.), O-ring P12 (2 pcs.), Stirring bar (2 pcs.)
Turntable connection cable (2m)	7075670K	For TTT-510 connection
AST connection cable (2m)	7075710K	For AST-3210 connection
RS-232 cable (2m)	0GC00002	For PC connection (D-sub9p for PC connection)
Reference electrode inner solution, 100mL (saturated KCI solution)	143F237	For reference electrode (HS-305DS and others)
Reference electrode outer solution RE-2, 100mL	143F238	For electrode for salinity titration
Printer chart paper (5 rolls/pack)	PAP-HCS	For internal printer, one of standard accessories of AUT-701

#### Parts for burette

Description	Code	Remarks
Unit plate (for 1mL to 20mL syringe)	7075610K	
Unit plate (for 50mL syringe)	7075620K	
20mL brown syringe	AUT-045P	
50mL brown syringe	P000010	
10mL brown syringe	AUT-046P	
5mL brown syringe	AUT-047P	
1mL brown syringe	AUT-048P	
Degassing nozzle (for 1mL to 20mL syringe)	P000070	
Degassing nozzle (for 50mL syringe)	P000071	
Titration nozzle (for micro titration)	AST-P008	
Reagent bottle with tube joint	7075630K	1,000mL translucent brown plastic bottle
Carbon dioxide adsorbent with tube	7075640K	One of standard accessories of ABT-7
Teflon tube, black, 2m (for 5mL to 20mL syringe)	AUT-022P	For 5mL to 50mL syringe
Teflon tube, black, 2m (for 1mL syringe)	AUT-024P	Only for 1mL syringe

#### Single Titration System Automatic Titrator AUT-701 + Automatic Burette ABT-7+ Magnetic Stirrer ST-7

#### **Automatic Titrator AUT-701**

Auto	matic	litrator AUI-/UI	
Display		Color graphic LCD with back light	
		Simultaneous indication of titration curve, titration volume, pH(mV),	
Displayed Item		temperature and titration time Changing over display of System A/System B display when dual	
		system	
		Ambient temperature and relative humidity when optional Temp. & RH sensor used	
Operatio	nal kev/	Flat kev	
	g method	Dialog method by function keys and ten keys	
		Inflection point detection titration	
		Set point titration	
		Full volume titration (inflection point detection/set point)	
T:44:	,	Intersection detection titration Stat titration	
Titration	ng method	Preliminary titration	
		pka measurement	
		pH/mV measurement	
		pH adjustment	
		Electric conductivity measurement	
Number of t	titration steps	Up to 5 steps (up to 2 steps for intersection titration)	
Titration	mada	Standard mode : 20/system User mode : 40/system	
Titration	mode	Link mode (for sequential titration): 20/system	
		Single titration system (system A)	
Titration	control	Dual titration system (system A/system B)	
system		Simultaneous titration/measurement available	
		pH : 0.00-14.00	
		mV : 0.0+/-2,000.0mV	
		Temperature : 0.0-100.0°C	
		Constant voltage current method : $0-2.000\mu$ A $0-20.00\mu$ A	
		0-20.00μA 0-200.0μA	
		0-1000 <i>μ</i> A	
		(by range selection)	
		Constant current voltage method : 0.0+/-2000.0mV	
		Electric conductivity : by cell  Displayed ranges are below: 0-200 0.45 /m(0-2 000 45 /m)	
Measurin	ng item/	Displayed ranges are below : $0-200.0\mu$ S/m( $0-2.000\mu$ S/m) $0-2.000$ mS/m( $0-20.00\mu$ S/m)	
range		0-2.000m3/m(0-20.00μ3/m) 0-20.00mS/m(0-200.0μS/m)	
		0-200.0mS/m(0-2.000mS/m)	
		0-2.000S/m(0-20.00mS/m)	
		0-20.00S/m(0-200.0mS/m)	
		0-200.0S/m(0-2.000S/m)	
		(Fixed range for titration, auto ranging for EC measurement)	
		Transmittance : 0.0-100.0%	
		(photometric titration)	
		Ambient temp. : 0.0-50.0℃ RH : 5-90%	
		Max. 4 channels	
Number		2 pcs. of pH/mV input installed as standard	
electrod	e inputs	3rd & 4th channels are for optional titration units	
	connectable	Max. 10 units	
auto. bure			
pH calibr	ation ard solution	Automatic 3-point calibration (manual calibration only for user's standard solutions)  JIS standard solutions/US standard solutions/User's specified standard solutions	
	and solution and sol setting	Given standard solutions table input: 3 points	
	compensation		
	measurement	0.0−100.0℃	
	Temperature	ATC (auto. temp. compensation) : 0.0-100.0°C	
Electric conductivity	compensation	MTC (manual temp. compensation): 0.0-100.0°C	
	range	TC off available	
compensation	Stand. temp. setting	0.0-100.0°C	
	Temp. coefficient	0.0-10.0%/°C	
		Ambient condition display, record (ambient temp., relative humidity) pH calibration history : Max. 20 calibrations	
		Self check up history : Max. 20 calibrations Self check up history : Max. 6 (by checking input)	
		Periodical check up history : Max. 10 (by standard solution)	
Validatio function	n support	Burette volume checking history: Max. 6 (data in burette)	
runction		Electrode management alarm (date management)	
		Reagent replacement alarm (date management)	
		Syringe replacement alarm (date or numbers of strokes management)	
Det-		Reagent empty alarm (alarm point setting available)	
Data me	mory is available)	100 data sets in case of single titration system 200 data sets in case of dual titration system	
Printer	aranabioj	Internal line thermal printer	
		RS-232 × 5ch.	
Communication		#1: Automatic burette (Max. 10 units connectable)	
		#2: PC	
		#3: Turntable (TTT-510)	
		#4: One from among Auto sampler AST-3210,	
		electric balance and Turntable TTT-510	
		#5: External printer EPS-R LAN × 1	
		USB(host) × 1	
Ambient temp. & RH		5-35°C, 85% or less (without condensation)	
Power so		AC100V-240V 50/60Hz	
Power consumption		Max. 35VA	
	ns & weight	Approx. W150 $\times$ H200 $\times$ D385 mm, Approx. 3.6 kg	

#### Standard Accessories

Description	Code
Magnetic stirrer	ST-7
pH electrode (Strong pH combined electrode)*	GST-5741C
pH4.01 standard solution, 500mL*	143F191
pH6.86 standard solution, 500mL*	143F192
3.3mol/L KCI solution, 50mL*	RE-4-05
Electrode stand with stopper	6948810K
Electrode holder	0IB00001
Electrode attachment (G)	0IB00004
Electrode attachment (J)	0IB00005
Electrode attachment (N)	0IB00008
Printer chart paper (2 rolls)	-
Instruction manual	-

<sup>\*</sup>Not attached when no electrode & no standard solution assigned

#### **Automatic Burette ABT-7**

Indication	O
iiiaioa tioii	Syringe size indication : LED
Syringe size setting	Digital rotary switch
Syringe	Fine inner grinding position type 20mL brown glass syringe 1mL, 5mL, 10mL and 50mL syringes are option
Wetted material	Teflon, Daiflon, hard glass
Using tube	ID2mm $\times$ OD3mm black Teflon tube (ID1mm $\times$ OD3mm black Teflon tube when 1mL syringe)
Absorption time	Approx. 20 seconds for full stroke
Discharge rate	50mL syringe : 2-150mL/min 20mL syringe : 0.6-60.0mL/min 10mL syringe : 0.3-30.0mL/min 5mL syringe : 0.2-15.0mL/min 1mL syringe : 0.03-3.00mL/min
Minimum dropping volume	50mL syringe: 0.0025mL 20mL syringe: 0.001mL 10mL syringe: 0.0005mL 5mL syringe: 0.00025mL 1mL syringe: 0.00005mL
Accuracy of burette	In case of 20mL syringe Whole volume error : +/-0.1% Repeatability : +/-0.01mL
Validation support	Burette volume checking history : Max. 6
Outlet for peripheral equipment	One for burette and one for magnetic stirrer
Ambient temp. & RH	5-35℃, 85% or less (without condensation)
Power source	AC100V-240V 50/60Hz
Power consumption	Approx. 30VA
Dimensions & weight	Approx. W110 $\times$ H396 $\times$ D346 mm, Approx. 5.4 kg

#### **Standard Accessories**

Description	Code
Unit plate (for 1mL to 20mL syringe)	7075610K
Reagent bottle with 1,000mL translucent brown plastic bottle	7075630K
Carbon dioxide adsorbent with tube	7075640K
20mL brown syringe	AUT-045P
Syringe attachment	AUT-066P
Degassing nozzle (for 1mL to 20mL syringe)	P000070
Nozzle cleaning attachment	AUT-067P
Nozzle purge tube	P000108
Teflon tube, black, 2m	AUT-022P
Burette connection cable, 0.9m	118B641
Instruction manual	-

## **Magnetic Stirrer**

Usable beaker	Up to 200mL beaker
External connection terminal	Power for overhead rod type stirrer
Power source	AC100V-240V 50/60Hz
Power consumption	Approx. 5VA
Dimensions & weight	Approx. W110 × H73 × D135 mm, Approx. 0.8 kg

#### **Standard Accessories**

Description	Code
Stirring bar	107D039

# Analysis Solutions for a Wide Range of Applications

# Chemistry



Composition analysis, saponification value, acid value, epoxy equivalence, diazotization titration

**Environmental** 

# Pharmaceuticals



Potentiometric titration, gastric acidity analysis

# **Food Processing**



Acidity, salinity, vitamin C, peroxide value

## Januar



Neutralization number and bromine number of petrochemical product, hydrogen sulfide in desulfurization waste water

# **Metal Plating**



Acid/alkali in nickel planting liquid, metal composition analysis

# Electronics

Acidity, alkalinity, hardness, residual chlorine



Free acid in etching liquid for printed circuit card, free acid in surface treatment liquid

# **DKK-TOA** CORPORATION



Do not operate producuts before consulting instruction manual.

International Operations: DKK-TOA Corporation

29-10, 1-Chome, Takadanobaba, Shinjuku-ku, Tokyo 169-8648 Japan Tel : +81-(0)3-3202-0225 Fax : +81-(0)3-3202-5685

Representative Office (Europe): DKK-TOA European Representative

St. Johns Innovation Centre, Cowley Rd., Cambridge CB4 0WS UK. Tel: +44 (0)1223-526471 Fax: +44 (0)1223-709239

http://www.dkktoa.net

Information and specifications are for a typical system and are subject to change without notice.